



INSTRUCTIONS FOR PATIENT USE



STEP 1:

Remove the protective cover with a light counterclockwise twist and pull while holding the barrel.



STEP 2:

Turn the dosing knob in the direction marked for the number of prescribed clicks.

NOTE: *Be sure to leave the knob in a "lock" position, and not between click positions.*



STEP 3:

Confirm that cream has exited the holes of the applicator cap.



STEP 4:

Apply the cream to the skin absorption site. **Rub** with the applicator cap until the cream is evenly spread (per physician or pharmacist direction). **Replace** the protective lid.



STEP 5:

Check the scale at the blue line for the number of clicks remaining. **Refill** your prescription as instructed by your pharmacy.

CONTACT YOUR PHARMACIST FOR MORE INFORMATION

EpitomeRx
(626) 631-6150
contact@epitomerx.com